



C₄ - Cut

#	Property	Units	Test Method	Value
1	1,3 Butadiene	wt%	G.C	43-60
2	Methyle Acetylene	wt%	G.C	0.15 max
3	Ethyl Acetylene	wt%	G.C	0.6 max
4	Vinyl Acetylene	wt%	G.C	2.5 max
5	1,2 Butadiene	wt%	G.C	0.8 max
6	I-Butane	wt%	G.C	0.3-4
7	N-Butane	wt%	G.C	2-9
8	Propadiene	wt%	G.C	0-1
9	T-2-Butene	wt%	G.C	2-7
10	1-Butene	wt%	G.C	8-13
11	I -Butene	wt%	G.C	15-26
12	C-2-Butene	wt%	G.C	1.5-5
13	C ₅ 's	wt%	G.C	<0.3

Shipping : Bulk Carrier Licensore: TPL - KTI